

Exergen TemporalScanner Calibration Choices



Exergen offers two options in using Temporal Artery Thermometers. The standard TemporalScanner is calibrated to provide a true core temperature, and the other is the same instrument, but recalibrated to display an oral equivalent temperature.

- **Standard:** The standard model, like a pulmonary artery catheter or esophageal probe, is measuring arterial (core) temperature. Arterial temperature is close to rectal temperature, approximately 0.8°F (0.4°C) higher than oral temperatures. This is by far the model selected by the vast majority of hospitals.
 - ✓ **An upward adjustment is required to account for the higher temperatures in order to avoid unnecessary fever workups (please see attached guidelines).**
- **Oral Equivalent:** The model with the oral equivalent calibration is programmed to compute the normal average cooling effect at the mouth, and automatically reduces the higher arterial temperature by that amount. This calibration allows the hospital to maintain existing protocols for fever workups based on oral temperature, and results in a reading consistent with the 98.6°F (37°C) mean normal oral temperature, in the range of 96.6 - 99.5°F (35.9 - 37.5°C) commonly seen with normal oral temperature measurements.
 - ✓ **No change in protocol is required if the hospital's protocol is based on an oral temperature.**

While oral temperature has never been considered a gold standard, and can easily mislead the clinician in identifying a fever, historically, it is the most common site for measuring body temperature. As such, the level of temperature requiring a fever/septic workup is typically based on an oral temperature, despite the fact that there can be a two degree difference (core/oral/axillary) among temperature measurement sites.

Nonetheless, it is not always easy to change collective memory (a normal temperature is still considered as 98.6°F (37°C) by the world in general (including clinicians), despite the site being measured, and despite the fact that only 8% of the world has a normal temp at 98.6°F (37°C). Accordingly, Exergen offers the oral equivalent calibration. The measurement is still as valid as the uncorrected standard model, but no change in protocol is required when using the oral equivalent model.

- Exergen TemporalScanner Thermometer, Model TAT-5000, Arterial (Standard) Model:
Part Number 124275
- Exergen TemporalScanner Thermometer, Model TAT-5000, Oral Equivalent Model:
Part Number 124375

For further information, please contact Exergen Customer Service:

Exergen Corporation, Watertown, MA

617-923-9900 x 6234

www.exergen.com -- www.TAThermometer.org

Adjustment in Fever Thresholds with Temporal Artery Temperature Assessment

Marybeth Pompei^a and Francesco Pompei, Ph.D.^{a,b}

Temporal artery temperature (TAT) is a core temperature, defined as the temperature of the blood perfusing the major organs. The physics and physiology of the measurement are designed to accomplish this by scanning the skin over the TA, then mathematically replacing the heat lost from the blood perfusion to the environment.¹ As a core temperature, TAT is confirmed to be comparable to pulmonary artery (PA) temperature,^{2,3,4} esophageal temperature,⁵ rectal temperature.^{6,7,8,9,10,11} As expected for a core temperature, TAT is approximately 0.4°C (0.8°F) higher than oral temperature.^{12,13} TAT has fewer and different local artifacts that cause variations than either oral or rectal temperatures, and thus will not always be exactly the same. In particular, oral temperature is heavily influenced by local cooling¹⁴, and rectal temperature has long been known to exhibit an inertia in response to thermal changes,^{15,16} recently observed in infants.¹⁷

The table below provides several fever guidelines and the recommended adjustment when using TAT. The adjustment rule is simply as follows:

- If current protocol is based on pulmonary artery, esophageal, or rectal temperatures: No change is necessary.
- If current protocol is based on oral temperatures: Add 0.4°C (0.8°F).

Fever Guidelines and Recommended Adjustment for Temporal Artery Temperature

Source	Fever Guidelines	Correction for Temporal Artery	Fever Guidelines Based on Temporal Artery
Infants and toddlers best practice guidelines ^{18,19,20,21,22}	1) infants 0 to 90 days of age with a RECTAL temperature $\geq 38.0^{\circ}\text{C}$ (100.4°F), and 2) infants and children 3 to 36 months of age with a RECTAL temperature $\geq 39.0^{\circ}\text{C}$ (102.2°F).	No change in existing protocol is necessary	1) infants 0 to 90 days of age with a RECTAL temperature $\geq 38.0^{\circ}\text{C}$ (100.4°F), and 2) infants and children 3 to 36 months of age with a RECTAL temperature $\geq 39.0^{\circ}\text{C}$ (102.2°F).
Infectious Diseases Society of America Guidelines – Revised 2002 ^{*23, 24}	A single ORAL temperature ≥ 38.3 (101°F) An ORAL temperature $\geq 38.0^{\circ}\text{C}$ (100.4°F) ≥ 1 hour	Add 0.4°C (0.8°F) Add 0.4°C (0.8°F)	A single TA temperature ≥ 38.7 (101.8°F) A TA temperature $\geq 38.4^{\circ}\text{C}$ (101.2°F) ≥ 1 hour
Principles and Practice of Pediatric Oncology ²⁵	One ORAL temperature $\geq 38.5^{\circ}\text{C}$ (101.3°F) Three ORAL temperatures $\geq 38.0^{\circ}\text{C}$ (100.4°F) at least 4 hours apart within 24 hours.	Add 0.4°C (0.8°F) Add 0.4°C (0.8°F)	One TA temperature ≥ 38.9 (102.1°F) Three TA temperatures $\geq 38.4^{\circ}\text{C}$ (101.2°F) at least 4 hours apart within 24 hours.
National Institutes of Health, Pediatric Oncology Branch, National Cancer Institute ^{**26}	One ORAL temperature $\geq 38.3^{\circ}\text{C}$ (101.0°F). Two oral temperatures $\geq 38.0^{\circ}\text{C}$ (100.4°F)	Add 0.4°C (0.8°F)	One TA temperature of $\geq 38.7^{\circ}\text{C}$ (101.8°F). Two TA temperatures $\geq 38.4^{\circ}\text{C}$ (101.2°F)

* Fever Definition in Neutropenic Patients with Unexplained Fever

** Evaluation and treatment of fever in the non-neutropenic child with cancer

^aExergen Corporation, 400 Pleasant Street, Watertown, MA, USA 02472; mpompei@exergen.com
^bDept. of Physics, Harvard University Cambridge, MA, USA 02139

-
- ¹ Pompei F, Pompei M. Non-invasive temporal artery thermometry: Physics, Physiology, and Clinical Accuracy, presented at *Medical Thermometry for SARS Detection, SPIE Defense and Security Symposium*, available in *Conference Proceedings*, April, 2004.
- ² Carroll DL, Finn C, Gill S, Sawyer J, Judge B. The Massachusetts General Hospital. A Comparison of measurements from a temporal artery thermometer and a pulmonary artery catheter thermistor. *AACN Poster* 2004.
- ³ Myny D, DeWaele J, Defloor T, Blot S, Colardyn F. Intensive Care Unit, Ghent University Hospital, Ghent, Belgium, Temporal scanner thermometry: a new method of core temperature measurement in intensive care patients. *SMJ* 2005 45(1): 15-18.
- ⁴ Lawson L, Bridges E, Ballou I, Eraker R, Greco S, Shively J, Sochulak V. University of Washington. Temperature measurement in critically ill adults. *Am. J. Crit. Care.*, May 2006; 15: 324 - 346.
- ⁵ Al-Mukazeem F, Allen U, Komar L, et al. University of Toronto/Hospital for Sick Children, Comparison of temporal artery, rectal and esophageal core temperatures in children: Results of a pilot study. *Journal of Pediatric and Child Health*, Vol 9, No 7, pp 461-465, 2004.
- ⁶ Greenes DS, Fleisher GR. Boston Childrens Hospital and Harvard Medical School. Accuracy of a noninvasive temporal artery thermometer for use in infants. *Arch Pediatr Adolesc Med*, Vol 155, pp 376-381, Mar 2001
- ⁷ Schuh S, Komar L, Stephens D, Chu L, Read S, Allen U, University of Toronto/Hospital for Sick Children. Comparison of the temporal artery and rectal thermometry in children in the emergency department. *Pediatric Academic Societies Annual Meeting*, May 3-6, 2003, Seattle, WA.
- ⁸ Schuh S, Komar L, Stephens D, Chu L, Read S, Allen U, University of Toronto/Hospital for Sick Children. Comparison of the temporal artery and rectal thermometry in children in the emergency department. *Pediatric Emergency Care*, Vol 20, No. 11, Nov 2004.
- ⁹ Siberry GK, Diener-West M, Schappell E, Karron RA, Department of Pediatrics, School of Medicine, The Johns Hopkins University. Comparison of temple temperatures with rectal temperatures in children under two years of age. *Clinical Pediatrics*, pp 405-414, July/August 2002.
- ¹⁰ Hebbbar K, Fortenberry JD, Rogers K, Merritt R, Easley K. Children's Healthcare of Atlanta at Egleston. Comparison of temporal artery thermometer to standard temperature measurements in pediatric intensive care unit patients. *Pediatr Crit Care Med*. 2005 Sep;6(5):557-61.
- ¹¹ Dybwik K, Nielsen EW. Infrared temporal temperature measurement. *Journal of the Norwegian Medical Association* 2003; 123: 3025-6.
- ¹² Routhier D, Hostler D, Wolfson A, Wheeler M, Reynolds J. University of Pittsburgh. Comparison of temporal artery and oral temperatures in the emergency department. *ACAD EMERG MED*, May 2006, Vol. 13, No. 5, Suppl. 1, www.aemj.org p. S99
- ¹³ Houdas Y, Ring, EFJ. *Human Body Temperature*. Ch 5, p83 Plenum Press, NY 1982.
- ¹⁴ Tandberg D, Sklar D. Effects of tachypnea on the estimation of body temperature by an oral thermometer. *NE J Med* v308(16):945-946m 1983.
- ¹⁵ Greenes DS, Fleisher GR. Boston Childrens Hospital and Harvard Medical School. When body temperature changes, does rectal temperature lag? *Journal of Pediatrics*, 02.037, pp 824-826, September 2004.
- ¹⁶ Houdas Y, Ring, EFJ. *Human Body Temperature*. Ch 5, p87 Plenum Press, NY 1982.
- ¹⁷ Cooper KE, Kenyon JR, A comparison of temperatures measured in the rectum, esophagus, and on the surface of the aorta during hypothermia in man, *Brit J Surgery*, pp 616-619, 1957.
- ¹⁸ Kourtis AP, Sullivan DT, and Sathian U, Practice Guidelines for the Management of Febrile Infants Less Than 90 Days of Age at the Ambulatory Network of a Large Pediatric Health Care System in the United States: Summary of New Evidence, *Clinical Pediatrics*, January 1, 2004; 43(1): 11 - 16.
- ¹⁹ Emergency Services: Algorithms to improve quality and patient satisfaction in the Emergency Department, National Association of Children's Hospitals and Related Institutions (NACHRI) Patient Care FOCUS Groups 2001.
- ²⁰ Baraff LJ, Schrager DL, Bass JW, Fleisher GR, Klein JO, McCracken GH Jr, Powell KR. Management of the young febrile child. Commentary on practice guidelines. *Pediatrics*. 1997 Jul;100(1):134-6.
- ²¹ Baraff LJ, Bass JW, Fleisher GR, Practice guideline for the management of infants and children 0 to 36 months of age with fever without source. *Pediatrics*. 1993; 92:1-12
- ²² Sharieff GQ, Hoecker C, Silva PD. *The Journal of Emergency Medicine*, Vol. 21, No. 1, pp. 1-6, 2001.
- ²³ Hughes WT, Armstrong D, Bodey GP, Brown AE, Edwards JE, Feld R, Pizzo P, Rolston KVI, Shenep JL, Young LS. 1997 Guidelines for the Use of Antimicrobial Agents in Neutropenic Patients with Unexplained Fever. Guidelines from the Infectious Diseases Society of America.
- ²⁴ Hughes WT, Armstrong D, Bodey GP, Bow EJ, Brown AE, Calandra T, Feld R, Pizzo PA, Rolston KV, Shenep JL, Young LS. 2002 guidelines for the use of antimicrobial agents in neutropenic patients with cancer. *Clin Infect Dis* 2002 Mar 15;34(6):730-51.
- ²⁵ Pizzo PA, Poplack DG. *Principles and Practice of Pediatric Oncology*, 2nd Ed, Part 6, Ch 42, Pg 1066. J.B. Lippincott Company, 1993.
- ²⁶ Salzer W, Steinberg SM, Liewehr DJ, Freifeld A, Balis FM, Widemann BC. Evaluation and Treatment of Fever in the Non-Neutropenic Child With Cancer. Pharmacology and Experimental Therapeutics Section, Pediatric Oncology Branch, Center for Cancer Research, National Cancer Institute, National Institutes of Health, *J Pediatr Hematol Oncol*. 2003 Aug; 25 (8):606-12.

